




PANOPTES COVID-19 presentation

Check !



Bring your forehead 30mm from point **A**.

Clear your forehead (hair is an obstacle) and remove your glasses, if you wear them.

The result is displayed in a few seconds on the screen (in **B**) and by a color in **C**.

Green : Everything is fine !

Red : Oops, you maybe got a cold ?
Consult government directives.

*The sensor measures the skin temperature, which varies between 33,5 and 36,9°C for a healthy subject.
The external temperature is lower than the internal body temperature, which varies between 36,1 and 37,8°C. The alarm is sent from 37°C.*

Thank you for participating

PANOPTES Team.

Visit the website www.panoptes-hm.com

Note :

PANOPTES is an indicative tool and doesn't replace medical advice. If you need more information and advice, ask your pharmacist. If symptoms worsen or persist, see a doctor.